


PulseCare

"Peace of Mind, Every Beat"



E/20/192 Karunaratne AGSI
E/20/418 Wahalathantri TN
E/20/440 Wickramsinghe RT



Problem

Cardia Care Gap



- Doctors only see patient data during hospital visits.
- Current care often relies on patients reporting symptoms after they occur, rather than capturing precautionary data to prevent incidents.
- Patient Isolation: High-risk patients often feel anxious and disconnected from medical oversight when recovering at home.

Solution



Core Concept: A compact wrist-worn device that continuously tracks Heart Rate and SpO₂ and bridges the gap between patient and guardian via BLE.


Key Value Proposition:

- **Continuous Monitoring:** Non-invasive tracking of vital signs.
- **Instant Alerts:** Emergency SOS button for immediate help.
- **Remote Connectivity:** Seamless data sync to a mobile app.

System Architecture

Input

MAX30100 Sensor
+
Button Input



Process

ESP32-C3 Super
Mini



Output

OLED Display
+
BLE Transmit



Receiver


Android App
+
Web Portal

android
studio 

 Cloud
Firestore

Demonstration

FireStore




Budget

Component	Unit Price (LKR)
3.7 V Li-Po ion battery	880
ESP32-C3 Super Mini	730
MAX30100 Sensor	290
TP4056 Charger IC	80
0.96" OLED 128×64 White	590
Button and Switch	60
Enclosure print	600
Package & Label design	300
Wires, Resistors, Capacitors, LED bulbs	100
Total Cost	3630

Heart Rate


- Uses Photoplethysmography (PPG) to detect blood volume changes in the microvascular bed of tissue.
- Normal Range: 60 – 100 BPM (Resting).
- Critical Detection: PulseCare identifies abnormal spikes (tachycardia) or drops (bradycardia) to trigger alerts.




SpO₂

SpO₂ measures the percentage of hemoglobin in the blood that is oxygenated. It acts as a direct report card for respiratory and cardiac efficiency.


- Healthy Range: 95% – 100%.
- Hypoxia Warning: Values below 90% indicate insufficient oxygen supply to the brain and organs.





White Box Analysis



Schematic



Enclosure Design



Thank you!

